
1

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
Reference
As observed by Mary Nabity, DVM, PhD, DAVCP,
Texas A&M University Publication pending

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Reference
Bratt J, Obare E, Yerranilli M, Elliott J, Yerranilli M. Relationship
between serum symmetric dimethylarginine concentration and
glomerular filtration rate in cats. J Vet Intern Med
2014;28(6):1699-1701.

5

International Renal Interest Society

[IRIS Guidelines](#) | [Education](#) | [Emerging Themes](#) | [About IRIS](#) | [Related sites](#)

Guidelines

- IRIS Staging of CKD
- IRIS Treatment Recommendations for CKD
- IRIS Grading of AKI
- IRIS Canine Proteinuria Consensus Project

IRIS Staging of CKD

As the leader in veterinary nephrology, IRIS has a responsibility to provide practicing veterinarians with evidence-based guidance for diagnosing, treating, and managing kidney disease. Our responsibility includes providing the most useful direction on how to interpret diagnostic information. The inclusion of SDMA values in the IRIS CKD Staging Guidelines is the next step in improving not only early diagnosis of CKD but also tailoring treatment of CKD to individual patients. The IRIS Board recommends using both serum creatinine and SDMA to enhance our ability to evaluate renal excretory function. These two surrogate markers are complementary to each other. Another change the IRIS board has recently made is the expansion of Stage 2 for dogs and a corresponding reduction in size of Stage 3 for dogs. This change was made to reduce the large variability in clinical signs and treatment needs for dogs in the previously large IRIS Stage 3 category.

Staging of chronic kidney disease (CKD) is undertaken following the diagnosis of CKD in order to facilitate appropriate treatment and monitoring of the patient. Staging is based initially on fasting blood creatinine, assessed on at least two occasions in the stable patient. The patient is then substage based on proteinuria and systemic blood pressure.

Based on these categories, some empirical recommendations can be made about the type of treatment it would be logical to use for these cases. In addition, predictions based on clinical experience might be made about the likely response to treatment.

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22

The Association Between Symmetric Dimethylarginine Concentrations and Neoplasia in Dogs and Cats

Michael Coyne, VMD PhD, Clinical Studies Design Specialist, IDEXX Laboratories, Inc.

The purpose of this study was to examine the association between SDMA concentrations in dogs and cats with different types of neoplasia. The IDEXX Laboratories pathology database was used to identify cases of lymphoma, hemangiosarcoma, mammary carcinoma, mammary adenocarcinoma, and lipoma in dogs, and lymphoma and visceral mast cell tumor in cats. Potential cases needed a serum SDMA and creatinine measured within 3 months prior to diagnosis. Control animals were randomly selected from the clinical laboratory database and age- and breed-matched (in dogs) and age-matched (in cats) to cases. SDMA and creatinine concentrations were compared between cases and controls by cancer type. SDMA concentrations were significantly higher in dogs and cats with lymphoma ($P < 0.0001$) and significantly lower in dogs with lipoma, mammary adenocarcinoma and mammary carcinoma ($P < 0.003$) as compared to controls. Creatinine concentrations were significantly lower in dogs with mammary adenocarcinoma and mammary carcinoma and in cats with lymphoma ($p < 0.0009$) when compared to controls. Odds ratios (OR) and 95% confidence intervals (CI) were calculated to determine the association between SDMA and tumor type (Table 1).

The results indicate that in both cats and dogs, cases of lymphoma are associated with increased SDMA concentrations with a large portion of these cases having concurrent creatinine concentrations within the reference interval. Although the etiology of increased SDMA concentrations in these cases was not determined,

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